

1.5T ADNI GE 12.0 M3 Software, TwinSpeed Gradient and 8-channel Brain Coil

Set up the patient according to the ADNI instruction manual. Use the 8-channel brain coil.

Use the electronic ADNI protocol that has been loaded on your scanner.

NOTE: The "adni" pulse sequence used for series 2 and 3 is a works in progress pulse sequence that ADNI sites obtain from GE Healthcare.

SERIES	1. 3 plane loc.	scan plane mode SAT	<i>IMAGING PARAMETERS</i>		<i>ACQUISITION TIMING</i>	
	coil 8hrbrain		3-plane (Whole Body gradient)	matrix/nex	256 / 128 / 1	
	etl			fov (cm)	26	
	scan time :13			slice/space autoshim	5 / 5 On	
comments	Use 8-channel brain coil.					

SERIES	2. Sag MPRAGE	scan plane mode pulse seq image opts. psd name	<i>IMAGING PARAMETERS</i>		<i>ACQUISITION TIMING</i>	
	coil 8hrbrain		Sag 3D (ZOOM gradient)	matrix/nex	192 / 192 / 1	
	#echos te		SPGR	phase fov		
	Prep time flip angle etl		EDR, IrPrep, Fast	locs/pause freq. direct.	S/l	
bw1/bw2 scan time	1000 8	adni	fc direct			
	<i>SCAN TIMING</i>	User CVs	<i>ADDITIONAL PARAMETERS</i>		phase corr	
	1		turbo mode = 1	autoshim	On	
	min full		swap fat chem shift = 1	<i>SCANNING RANGE</i>		
	15.63		slice resolution = 100%	fov	24	
	7:42	mprage mode = 1	slice/space	1.2mm 170 locs/slab		
		TD or mprageTR=2400				
comments	Cover skin to skin. Remind the patient to hold still for this scan.					

SERIES	3.MPRAGE-repeat	scan plane mode pulse seq image opts. psd name	<i>IMAGING PARAMETERS</i>		<i>ACQUISITION TIMING</i>	
	coil 8hrbrain		Sag 3D (ZOOM gradient)	matrix/nex	192 / 192 / 1	
	#echos te		SPGR	phase fov		
	Prep time flip angle etl		EDR, IrPrep, Fast	locs/pause freq. direct.	S/l	
bw1/bw2 scan time	1000 8	adni	fc direct			
	<i>SCAN TIMING</i>	User CVs	<i>ADDITIONAL PARAMETERS</i>		phase corr	
	1		turbo mode = 1	autoshim	Off	
	min full		swap fat chem shift = 1	<i>SCANNING RANGE</i>		
	15.63		slice resolution = 100%	fov	24	
	7:42	mprage mode = 1	slice/space	1.2mm 170 locs/slab		
		TD or mprageTR=2400				
comments	Prescribe same image locations as series 2, unless adjustment is needed (e.g. to correct for wrap). Remind the patient to hold still for this scan.					

(Continued on next page)

1.5T ADNI GE 12.0 M3 Software, TwinSpeed Gradient and 8-channel Brain Coil

SERIES coil	4. Sag B1 Cal PA	scan plane mode	Sag	matrix/nex	128 / 128 / 1
	8hrbrain		3D (ZOOM gradient)		phase fov
	<i>SCAN TIMING</i>	pulse seq	GRE	locs/pause	
#echos	1	image opts.	EDR, Fast	freq. direct.	S/l
te	min full	psd name		fc direct	
tr			<i>ADDITIONAL PARAMETERS</i>	Auto Shim	On
flip angle	2				<i>SCANNING RANGE</i>
etl		User CVs	Image acq. delay=0	fov	30
bw1/bw2	62.5		Turbo Mode =1	slice/space	2.5mm locs/slab=96
scan time	:36		Slice resolution=100%		

comments Cover skin to skin.

SERIES coil	5. Sag B1 Cal PA	scan plane mode	Sag	matrix/nex	128 / 128 / 1
	BODY		3D (ZOOM gradient)		phase fov
	<i>SCAN TIMING</i>	pulse seq	GRE	locs/pause	
#echos	1	image opts.	EDR, Fast	freq. direct.	S/l
te	min full	psd name		fc direct	
tr			<i>ADDITIONAL PARAMETERS</i>	Auto Shim	Off
flip angle	2				<i>SCANNING RANGE</i>
etl		User CVs	Image acq. delay=0	fov	30
bw1/bw2	62.5		Turbo Mode =1	slice/space	2.5mm locs/slab=96
scan time	:36		Slice resolution=100%		

comments Turn Body Coil on for this scan. Leave 8hrbrain coil plugged in, but accept change to Body Coil
Same image locations as series 4.

SERIES coil	6. Ax PD/T2 FSE	scan plane mode	Ax	matrix/nex	256 / 256 / 1
	8hrbrain		2D (ZOOM gradient)		phase fov
	<i>SCAN TIMING</i>	pulse seq	FSE-XL	acqs/pause	0
#echoes	2	image opts.	EDR, Fast	freq. direct.	A/P
te	min full / TE2=100	psd name		fc direct	
TR	3000		<i>ADDITIONAL PARAMETERS</i>	Autoshim	Off
flip angle				phase corr	
etl	16	User CVs	Filter Choices=none		<i>SCANNING RANGE</i>
bw1/bw2	20.83		blurring cancellation=0	fov	24
scan time	6:25			slice/space	48 loc, 3mm interleaved

comments Accept change back to 8-channel Brain Coil
Prescribe 48 slices to cover head.

Series 6 is the final patient series.

Follow the ADNI instructions to complete the phantom scans.