

**1.5T ADNI GE 12.0 M4 Software, TwinSpeed Gradient and 8-channel Brain Coil**

Set up the patient according to the ADNI instruction manual. Use the 8-channel brain coil.

Use the electronic ADNI protocol that has been loaded on your scanner.

**NOTE: The "adni" pulse sequence used for series 2 and 3 is a works in progress pulse sequence that ADNI sites obtain from GE Healthcare.**

<b>SERIES</b>	1. 3 plane loc.	scan plane mode SAT	<i>IMAGING PARAMETERS</i>		<i>ACQUISITION TIMING</i>	
	coil 8hrbrain		3-plane (Whole Body gradient)	matrix/nex	256 / 128 / 1	
	etl			fov (cm)	26	
	scan time		:13	slice/space	5 / 5	
			autoshim	On		
comments	Use 8-channel brain coil.					

<b>SERIES</b>	2. Sag MPRAGE	scan plane mode pulse seq image opts. psd name	<i>IMAGING PARAMETERS</i>		<i>ACQUISITION TIMING</i>	
	coil 8hrbrain		Sag 3D ( ZOOM gradient)	matrix/nex	192 / 192 / 1	
	#echos		1	phase fov		
	te		min full	locs/pause		
Prep time	1000	EDR, IrPrep, Fast	freq. direct.	S/I		
flip angle	8	adni	fc direct			
etl		<i>ADDITIONAL PARAMETERS</i>	phase corr			
bw1/bw2	15.63	Image acq. delay = 0	autoshim	On		
scan time	7:42	turbo mode = 1				
		swap fat chem shift = 1	fov	24		
		slice resolution = 100%	slice/space	1.2mm 170 locs/slab		
		mprage mode = 1				
		TD or mprageTR=2400				
comments	Cover skin to skin. Remind the patient to hold still for this scan.					

<b>SERIES</b>	3.MPRAGE-repeat	scan plane mode pulse seq image opts. psd name	<i>IMAGING PARAMETERS</i>		<i>ACQUISITION TIMING</i>	
	coil 8hrbrain		Sag 3D (ZOOM gradient)	matrix/nex	192 / 192 / 1	
	#echos		1	phase fov		
	te		min full	locs/pause		
Prep time	1000	EDR, IrPrep, Fast	freq. direct.	S/I		
flip angle	8	adni	fc direct			
etl		<i>ADDITIONAL PARAMETERS</i>	phase corr			
bw1/bw2	15.63	Image acq. delay = 0	autoshim	Off		
scan time	7:42	turbo mode = 1				
		swap fat chem shift = 1	fov	24		
		slice resolution = 100%	slice/space	1.2mm 170 locs/slab		
		mprage mode = 1				
		TD or mprageTR=2400				
comments	Prescribe same image locations as series 2, unless adjustment is needed (e.g. to correct for wrap). Remind the patient to hold still for this scan.					

**NOTE: Be sure to select the "Copy FOV, Thickness, Spacing" button for the graphic prescription. Otherwise, if you copy the slice locations from series 2, the number of slices might be reduced.**

(Continued on next page)

## 1.5T ADNI GE 12.0 M4 Software, TwinSpeed Gradient and 8-channel Brain Coil

<b>SERIES</b> coil	4. Sag B1 Cal PA	scan plane mode	Sag	matrix/nex	128 / 128 / 1
	8hrbrain		3D (ZOOM gradient)		phase fov
	<i>SCAN TIMING</i>	pulse seq	GRE	locs/pause	
#echos	1	image opts.	EDR, Fast	freq. direct.	S/I
te	min full	psd name		fc direct	
tr			<i>ADDITIONAL PARAMETERS</i>	Auto Shim	On
flip angle	2				<i>SCANNING RANGE</i>
etl		User CVs	Image acq. delay = 0	fov	30
bw1/bw2	62.5		Turbo Mode = 1	slice/space	2.5mm locs/slab=96
scan time	:36		Slice resolution = 100%		

comments Cover skin to skin.

<b>SERIES</b> coil	5. Sag B1 Cal PA	scan plane mode	Sag	matrix/nex	128 / 128 / 1
	<b>BODY</b>		3D (ZOOM gradient)		phase fov
	<i>SCAN TIMING</i>	pulse seq	GRE	locs/pause	
#echos	1	image opts.	EDR, Fast	freq. direct.	S/I
te	min full	psd name		fc direct	
tr			<i>ADDITIONAL PARAMETERS</i>	Auto Shim	Off
flip angle	2				<i>SCANNING RANGE</i>
etl		User CVs	Image acq. delay=0	fov	30
bw1/bw2	62.5		Turbo Mode =1	slice/space	2.5mm locs/slab=96
scan time	:36		Slice resolution = 100%		

comments Turn Body Coil on for this scan. Leave 8hrbrain coil plugged in, but accept change to Body Coil (In other words, select "**Apply**".) Same image locations as series 4.

<b>SERIES</b> coil	6. Ax PD/T2 FSE	scan plane mode	Ax	matrix/nex	256 / 256 / 1
	8hrbrain		2D (ZOOM gradient)		phase fov
	<i>SCAN TIMING</i>	pulse seq	FSE-XL	acqs/pause	0
#echoes	2	image opts.	EDR, Fast	freq. direct.	A/P
te	min full / TE2=100	psd name		fc direct	
TR	3000		<i>ADDITIONAL PARAMETERS</i>	Autoshim	Off
flip angle				phase corr	
etl	16	User CVs	blurring cancellation=0		<i>SCANNING RANGE</i>
bw1/bw2	20.83			fov	24
scan time	6:25			slice/space	48 loc, 3mm interleaved

comments Accept change back to 8-channel Brain Coil. (In other words, select "**Apply**".) Prescribe 48 slices to cover head.

**Series 6 is the final patient series.**

**Follow the ADNI instructions to complete the phantom scans.**