

1.5T ADNI GE 9.1 M4 Software, QC Phantom
BRM Gradient and a Single-Channel Birdcage Head Coil

Set up the phantom according to the ADNI instruction manual. Use the transmit/receive head coil

Use the electronic ADNI QC Phantom protocol that has been loaded on your scanner.

NOTE: The "adni" pulse sequence used for series 2 is a works in progress pulse sequence that ADNI sites obtain from GE Healthcare.

SERIES coil etl scan time	1. 3 plane loc.	scan plane mode SAT	<i>IMAGING PARAMETERS</i>		<i>ACQUISITION TIMING</i>	
	head		3-plane	matrix/nex	256 / 128 / 1	
				fov (cm)	26	
	:19			slice/space	5 / 5	
comments	Use the birdcage head coil		autoshim	On		

SERIES coil #echos te Prep time flip angle etl bw1/bw2 scan time	2. Sag QC MPRAGE	scan plane mode pulse seq image opts. psd name	<i>IMAGING PARAMETERS</i>		<i>ACQUISITION TIMING</i>	
	head		Sag	matrix/nex	192 / 192 / 1	
			3D	phase fov		
			SPGR	locs/pause		
		SCAN TIMING	EDR, IrPrep, Fast	freq. direct.	S/I	
	1		adni	fc direct		
	min full		<i>ADDITIONAL PARAMETERS</i>	phase corr		
	1000	User CVs	swap fat chem shift = 1	autoshim	On	
	8		special = 0			
			burst = 0	fov	24	
15.63		mprage mode =1	slice/space	1.3mm 184 locs/slab		
9:36		TD or mprageTR=3000		<i>SCANNING RANGE</i>		
comments	Cover the phantom.					