

SERIES 1. 3 Plane Localizer

Patient (RF)		Imaging		Acquisition Timing	
Coil	8HRBRAIN	Scan plane	3-PLANE	Frequency	256
Scan Timing		Mode	2D	Phase	128
		Pulse sequence	Gradient Echo	NEX	1.00
		Imaging Options	Seq, Fast	phase field of view	1.00
Additional Parameters				phase correct	No
				swap phase & frequency	Unswap
				autoshim	Auto
Range / Prescription					
				field of view	26.0
				slice thickness	5.0
				slice spacing	0.0

SERIES 2. Sag IR-FSPGR

Patient (RF)		Imaging		Acquisition Timing	
Coil	8HRBRAIN	Scan plane	SAGITTAL	Frequency	256
Scan Timing		Mode	3D	Phase	256
TE	Min Full	Pulse sequence	SPGR	NEX	1.00
TI	400	Imaging Options	EDR, Fast, IrP	phase field of view	1.00
Flip angle	11	Additional Parameters			
BW1	31.25	User CVs		phase correct	No
		CV4	0.00	swap phase & frequency	S/I
		CV5	0.00	autoshim	On
		CV6	1.00	Range / Prescription	
		CV23	100.00	slices per slab	200
				field of view	26.0
				slice thickness	1.2
				number of slices	1

Comments Please do not oblique. Please do not reduce the number of slices.

SERIES 3. Calibration Scan

Patient (RF)		Imaging		Acquisition Timing	
Coil	8HRBRAIN	Scan plane	AXIAL	phase correct	No
Scan Timing		Mode	2D	swap phase & frequency	R/L
		Pulse sequence	Gradient Echo	autoshim	Auto
		Imaging Options	Fast, Calib	Range / Prescription	
Additional Parameters				field of view	30.0
				slice thickness	5.0
				slice spacing	0.0
				number of slices	60

Comments Please do not oblique. Please do not reduce the number of slices.

SERIES 4. Accelerated Sag IR-FSPGR

Patient (RF)		Imaging		Acquisition Timing	
Coil	8HRBRAIN	Scan plane	SAGITTAL	Frequency	256
Scan Timing		Mode	3D	Phase	256
TE	Min Full	Pulse sequence	SPGR	NEX	1.00
TI	400	Imaging Options	EDR, Fast, IrP, Asset 	phase field of view	1.00
Flip angle	11	Additional Parameters		phase correct	No
BW1	31.25	User CVs		swap phase & frequency	S/I
		CV4	0.00	autoshim	On
		CV5	0.00	Range / Prescription	
		CV6	1.00	slices per slab	200
		CV23	100.00 	field of view	27.0
				slice thickness	1.2
				number of slices	1

Comments Please do not oblique. Please do not reduce the number of slices.

SERIES 5. Axial FLAIR

Patient (RF)		Imaging		Acquisition Timing	
Coil	8HRBRAIN	Scan plane	AXIAL	Frequency	256
Scan Timing		Mode	2D	Phase	256
TE	147.0	Pulse sequence	T2flair	NEX	1.00
TR	11000.0	Imaging Options	Fast	phase correct	No
TI	2250	Additional Parameters		swap phase & frequency	A/P
BW1	31.25	User CVs		autoshim	Auto
		CV3	3.00	Range / Prescription	
		CV8	0.00	field of view	22.0
		CV9	0.00	slice thickness	5.0
		CV11	0.00	slice spacing	0.0
		CV14	0.00	number of slices	42
		CV21	0.00		
		CV22	0.00 		

Comments Please do not oblique. Please do not reduce the number of slices. 

SERIES 6. Axial T2 Star

Patient (RF)		Imaging		Acquisition Timing	
Coil	8HRBRAIN	Scan plane	AXIAL	Frequency	256
Scan Timing		Mode	2D	Phase	192
Number of echoes	1	Pulse sequence	Gradient Echo	NEX	1.00
TE	20.0	Imaging Options	FC, EDR	phase field of view	1.00
TR	650.0	Additional Parameters		phase correct	No
Flip angle	20			swap phase & frequency	A/P
BW1	25.00			autoshim	Auto
				Range / Prescription	
				field of view	20.0
				slice thickness	4.0
				slice spacing	0.0
				number of slices	50

Comments Please do not oblique. Please do not reduce the number of slices.

SERIES 7. Axial DTI

Patient (RF)		Imaging		Acquisition Timing	
Coil	8HRBRAIN	Scan plane	AXIAL	Frequency	128
Scan Timing		Mode	2D	Phase	128
TE	Minimum	Pulse sequence	Spin Echo	phase field of view	0.66
TR	9050.0	Imaging Options	EPI, DIFF, Asset 	phase correct	Yes
		Additional Parameters		swap phase & frequency	R/L
				autoshim	Auto
		User CVs		Range / Prescription	
		CV5	1.00	field of view	35.0
		CV9	0.00	slice thickness	2.7
		CV18	0.00	slice spacing	0.0
				number of slices	59

Comments Please do not oblique. Please do not reduce the number of slices.

If you notice that the slice thickness is changed from 2.7 to 2.9 you can fix this bug by turning fat sat off then turning it back on again.



SERIES 8. Ax T2 FSE with Fat Sat**Patient (RF)**

Coil 8HRBRAIN

Scan Timing

Number of echoes 1
 TE 80.0
 TR 4000.0
 ETL 12
 BW1 19.23

Imaging

Scan plane AXIAL
 Mode 2D
 Pulse sequence FRFSE-XL
 Imaging Options FC, TRF, Fast, FR

Additional Parameters**User CVs**

CV7 0.00
 CV8 0.00
 CV14 0.00
 CV19 1.00
 CV21 1.00
 CV22 1.00

Acquisition Timing

Frequency 256
 Phase 256
 NEX 1.00
 phase field of view 1.00
 phase correct Yes
 Flow comp direction Slice
 swap phase & frequency A/P
 autoshim Auto

Range / Prescription

field of view 24.0
 slice thickness 4.0
 slice spacing Intleave
 number of slices 44

Comments Please do not oblique. Please do not reduce the number of slices.