

```

Nucleus = "H1";
SmartSelect = "yes";
Coil 1 (exclude) = "None";
Uniformity = "CLEAR";
FOV      FH (mm) = 256;
         AP (mm) = 256;
         RL (mm) = 192;
Voxel size  FH (mm) = 1;
            AP (mm) = 1.00392163;
            RL (mm) = 1.20000005;
Recon voxel size (mm) = 1;
Fold-over suppression = "no";
Slice oversampling = "user defined";
oversample factor = 1;
RF select. FOS = "no";
Reconstruction matrix = 256;
SENSE = "yes";
P reduction (AP) = 3;
S reduction (RL) = 1;
k-t BLAST = "no";
Overcontiguous slices = "no";
Stacks = 1;
slices = 160;
slice orientation = "sagittal";
fold-over direction = "AP";
fat shift direction = "F";
Stack Offc. AP (P=+mm) = 6.04466534;
            RL (L=+mm) = 0;
            FH (H=+mm) = -2.67916512;
Ang. AP (deg) = 0;
     RL (deg) = -0;
     FH (deg) = -0;
Free rotatable = "no";
Multi-chunk = "no";
Large table movement = "no";
PlanAlign = "no";
REST slabs = 0;
Catheter tracking = "no";
Interactive positioning = "no";
Allow table movement = "no";
Patient position = "head first";
orientation = "supine";
Scan type = "Imaging";
Scan mode = "3D";
technique = "IR";
Acquisition mode = "cartesian";
Fast Imaging mode = "TSE";
3D VIEW = "Brain FLAIR";
shot mode = "multishot";
TSE factor = 150;
startup echoes = 6;
profile order = "linear";
turbo direction = "Y";
DRIVE = "no";
fid reduction = "default";
Echoes = 1;
partial echo = "no";
TE = "user defined";
(ms) = 271;
Refocusing control = "yes";
angle (deg) = 40;
echo enhancement = "no";
bright fat reduction = "no";

```

```

TR = "user defined";
  (ms) = 4800;
Halfscan = "no";
Water-fat shift = "user defined";
  (pixels) = 0.519999981;
IR delay (ms) = 1650;
  dual = "no";
  power = "1";
RF Shims = "fixed";
Shim = "default";
mDIXON = "no";
Fat suppression = "SPIR";
  strength = "strong";
  frequency offset = "default";
Grad Rev Fat suppression = "no";
Water suppression = "no";
MTC = "no";
T2prep = "yes";
  echo time (ms) = 125;
  refocusing pulses = 4;
  flow venc (cm/s) = (3) 0;
Research prepulse = "no";
Zoom imaging = "no";
Diffusion mode = "no";
Elastography mode = "no";
Transmit channels = "both";
SAR mode = "high";
B1 mode = "default";
SAR Patient data = "auto";
PNS mode = "high";
Gradient mode = "maximum";
SoftTone mode = "no";
Cardiac synchronization = "no";
Heart rate > 250 bpm = "no";
Respiratory compensation = "no";
Navigator respiratory comp = "no";
Flow compensation = "no";
Motion smoothing = "no";
NSA = 1;
CENTRA = "no";
Manual start = "no";
Dynamic study = "no";
Arterial Spin labeling = "no";
Preparation phases = "auto";
Interactive F0 = "no";
SmartPlan survey = "no";
B0 field map = "no";
B1 field map = "no";
MIP/MPR = "no";
Images = "M", (3) "no";
Autoview image = "M";
Reference tissue = "Grey matter";
Recon compression = "No";
Preset window contrast = "soft";
Reconstruction mode = "real time";
Save raw data = "no";
Hardcopy protocol = "no";
Image filter = "system default";
Geometry correction = "default";
Elliptical k-space shutter = "default";
IF_info_seperator = 1634755923;
Total scan duration = "05:45.6";
Rel. SNR = 1.10176492;

```

Act. TR/TI (ms) = "4800 / 1650";  
Act. TE (ms) = "271";  
ACQ matrix M x P = "256 x 255";  
ACQ voxel MPS (mm) = "1.00 / 1.00 / 1.20";  
REC voxel MPS (mm) = "1.00 / 1.00 / 1.20";  
Scan percentage (%) = 100;  
WFS (pix) / BW (Hz) = "0.454 / 957.4";  
TSE es / shot (ms) = "3.2 / 513";  
TEeff / TEequiv (ms) = "271 / 119";  
Min. TR/TI (ms) = "2312 / 50";  
SAR / local torso = "< 15 %";  
Whole body / level = "< 0.5 W/kg / normal";  
SED = "< 0.2 kJ/kg";  
B1+rms / Coil Power = "0.88 uT / 14 %";  
Max B1+rms = "0.97 uT";  
PNS / level = "81 % / 1st level";  
dB/dt = "108.5 T/s";  
Sound Pressure Level (dB) = 27.2498283;  
Boil-off (hPa/h) = 1.63544679;  
Max Boil-off (hPa/h) = 15.279273;  
Boil-off (hPa) = 0.157002896;

```

Nucleus = "H1";
SmartSelect = "yes";
Coil 1 (exclude) = "None";
Uniformity = "CLEAR";
FOV RL (mm) = 256;
AP (mm) = 256;
FH (mm) = 160;
Voxel size RL (mm) = 2;
AP (mm) = 2;
Slice thickness (mm) = 2;
Recon voxel size (mm) = 2;
Small FOV imaging = "no";
Fold-over suppression = "no";
Reconstruction matrix = 128;
SENSE = "yes";
P reduction (AP) = 2;
k-t BLAST = "no";
Stacks = 1;
type = "parallel";
slices = 80;
slice gap = "user defined";
gap (mm) = 0;
slice orientation = "transverse";
fold-over direction = "AP";
fat shift direction = "P";
Stack Offc. AP (P=+mm) = 6.04466534;
RL (L=+mm) = 0;
FH (H=+mm) = -2.67916512;
Ang. AP (deg) = 0;
RL (deg) = -0;
FH (deg) = -0;
Free rotatable = "no";
Minimum number of packages = 1;
Slice scan order = "default";
Large table movement = "no";
PlanAlign = "no";
REST slabs = 0;
Catheter tracking = "no";
Interactive positioning = "no";
Allow table movement = "no";
Patient position = "head first";
orientation = "supine";
Scan type = "Imaging";
Scan mode = "MS";
technique = "SE";
Modified SE = "no";
Acquisition mode = "cartesian";
Fast Imaging mode = "EPI";
shot mode = "single-shot";
Echoes = 1;
partial echo = "no";
TE = "shortest";
Flip angle (deg) = 90;
TR = "shortest";
Halfscan = "yes";
factor = 0.868852437;
Water-fat shift = "minimum";
RF Shims = "fixed";
Shim = "auto";
mDIXON = "no";
Fat suppression = "SPIR";
strength = "strong";
frequency offset = "user defined";

```

```

offset (Hz) = 175;
Grad Rev Fat suppression = "no";
Water suppression = "no";
BB pulse = "no";
MTC = "no";
Research prepulse = "no";
Diffusion mode = "DTI";
    sequence = "SE";
    gradient duration = "maximum";
    gradient overplus = "no";
    directional resolution = "from file";
    average high b = "no";
Elastography mode = "no";
Transmit channels = "both";
SAR mode = "high";
B1 mode = "default";
SAR Patient data = "auto";
PNS mode = "moderate";
Gradient mode = "maximum";
SoftTone mode = "no";
Cardiac synchronization = "no";
Heart rate > 250 bpm = "no";
Respiratory compensation = "no";
Navigator respiratory comp = "no";
Flow compensation = "no";
Temporal slice spacing = "default";
NSA = 1;
Manual start = "no";
Dynamic study = "no";
    dyn stabilization = "regular";
Arterial Spin labeling = "no";
Preparation phases = "auto";
Interactive F0 = "no";
SmartPlan survey = "no";
B0 field map = "no";
B1 field map = "no";
MIP/MPR = "no";
Images = "M", (3) "no";
Autoview image = "M";
Calculated images = (4) "no";
Reference tissue = "White matter";
EPI 2D phase correction = "no";
Recon compression = "No";
Preset window contrast = "soft";
Reconstruction mode = "immediate";
Save raw data = "no";
Hardcopy protocol = "no";
Image filter = "system default";
Geometry correction = "default";
IF_info_seperator = 1634755923;
Total scan duration = "08:05.1";
Rel. SNR = 1.00341499;
Act. TR (ms) = "10861";
Act. TE (ms) = "100";
ACQ matrix M x P = "128 x 126";
ACQ voxel MPS (mm) = "2.00 / 2.03 / 2.00";
REC voxel MPS (mm) = "2.00 / 2.00 / 2.00";
Scan percentage (%) = 98.4375;
Packages = 1;
Min. slice gap (mm) = -0;
Diffusion gradient timing DELTA / delta (ms) = "50.0 / 25.9";
User defined DTI scheme =

```

```
"ADNI3_Basic_Philips7 (36, 100";
EPI factor = 63;
WFS (pix) / BW (Hz) = "19.393 / 22.4";
BW in EPI freq. dir. (Hz) = "1865.9";
SPIR offset act./default (Hz) = "175 [220]";
SAR / local torso = "< 30 %";
Whole body / level = "< 1.0 W/kg / normal";
SED = "< 0.5 kJ/kg";
B1+rms / Coil Power = "1.24 uT / 28 %";
Max B1+rms = "1.24 uT";
PNS / level = "79 % / normal";
dB/dt = "88.0 T/s";
Sound Pressure Level (dB) = 19.0595226;
Boil-off (hPa/h) = 12.4823666;
Max Boil-off (hPa/h) = 12.4826345;
Boil-off (hPa) = 1.68187499;
```