3.0T ADNI-Related GE 14.0 M4 Software, TwinSpeed Gradient and 8-channel Brain Coil

(Not for patient use--this protocol is for use with the ADNI phantom)

This protocol provides suggested imaging parameters for research studies that want to approximate the imaging methods used in the ADNI study using 14.0 M4 HDx software at 3T, but do not have access to an MP-RAGE pulse sequence. For more details see the document: "Use of ADNI MRI Methods for Non-ADNI Studies".

Accept the "First Operating Mode" pop-up in Series 1. Consult the scanner's user's manual to understand this choice and its implications.

*The head portion of the 16-channel head-neck-spine (HNS) can be used instead of the 8-channel brain coil for the entire study, if desired.

			IMAGING PARAMETERS		ACQUISITION TIMING
SERIES	1. 3 plane loc.	scan plane	3-plane	matrix/nex	256 / 128 / 1
coil	8hrbrain*	mode	(Whole Body gradient)	fov (cm)	26
etl		SAT		slice/space	5/5
scan time	:13			autoshim	On

comments	Use 8-channel brain coil*				
SERIES	2. Calibration Scan.	scan plane	axial	matrix/nex	default
coil	8hrbrain*	mode	(Whole Body gradient)	fov (cm)	30
etl		SAT	, , , , , , , , , , , , , , , , , , , ,	slice/space	6/0 43 slices
scan time	:13			autoshim	off
		•••			
comments	Used for PURE B1-corre	ction. Be sure t	o cover phantom completely	′.	
	(Although 14.0 M4 does i	not support PU	RE for 3T, this acquisition of	the PURE cali	bration scan
	is recommended for retro	spective proce	ssing in case PURE become	es compatible in	n future releases.)
			IMAGING PARAMETERS		ACQUISITION TIMING
SERIES	3. Sag IR-FSPGR	scan plane	Sag	matrix/nex	256 / 256 / 1
coil	8hrbrain*	mode	3D (ZOOM gradient)	phase fov	0.94
	SCAN TIMING	pulse seq	SPGR	locs/pause	
#echos	1	image opts.	EDR, IrPrep, Fast	freq. direct.	S/I
te	min full	psd name	efgre3d_cs	fc direct	
Prep time	650		ADDITIONAL PARAMETERS	phase corr	
flip angle	8		Image acq. delay = 0	autoshim	On
etl		User CVs	Turbo mode = 1		SCANNING RANGE
bw1/bw2	31.25		Slice resolution = 100%	fov	26
scan time	9:43			slice/space	1.3mm 170 locs/slab
	h	ool	<u> </u>	·	
					L
comments	Cover phantom complete	ely			