

Notice to site personnel uploading ADNI MR and PET scans

A system modification will be implemented on Thursday, April 4th 2013. When uploading subject scan data, you will no longer need to select a research group from the drop-down listing (Patient, Phantom, Volunteer). Instead, only phantom and volunteer uploads will require selection from a drop-down list.

Before change

Please follow the instructions outlined above:

Project ADNI@BostonU Bypass validation steps

Select Data Type Original XML

Research Group

Visit

Subject ID: Max. 10 characters allowed
Identifier to replace Patient ID

Source Directory: BROWSE...
Location of original files

Target Directory: BROWSE...
Location for target files

NOTE: Source Directory for file formats with complete headers (DICOM, GE, ECAT, etc) may contain multiple series from a single subject. Source Directory for file formats with limited headers (ANALYZE, MINC) or no headers (TIFF, TGA, etc) must contain a single series for a single subject in a single directory. Headerless files must contain a sequential slice number within the file name. ANALYZE files are assumed to be in SPM orientation.

Record diagnostics to file

After change

Please follow the instructions outlined above:

Project ADNI@BostonU Bypass validation steps

Select Data Type Original XML

Visit

Subject ID: Max. 10 characters allowed
Identifier to replace Patient ID

Source Directory: BROWSE...
Location of original files

Target Directory: BROWSE...
Location for target files

NOTE: Source Directory for file formats with complete headers (DICOM, GE, ECAT, etc) may contain multiple series from a single subject. Source Directory for file formats with limited headers (ANALYZE, MINC) or no headers (TIFF, TGA, etc) must contain a single series for a single subject in a single directory. Headerless files must contain a sequential slice number within the file name. ANALYZE files are assumed to be in SPM orientation.

Record diagnostics to file

For non-subject scans, a pop-up window will appear to offer research group selection:

Select a Research Group

Please select the research group